


## BURTONS LIFETIME STAINLESS STEEL CAGE BANK FOR DOG/MIXED WARD

170-028-05



### PRODUCT DESCRIPTION

## Built from the highest quality materials, Burtons Lifetime Cages are guaranteed to last a lifetime!

This pre-configured cage bank consists of 5 cages, mesh doors and required fixings to make up a bank included & is designed to be used for Dog/mixed ward.

#### Includes:

**1 x 36"H x 48"W double door stainless steel lifetime cage**

**1 x 36"H x 36"W stainless steel lifetime cage**

**1 x 30"H x 24"W stainless steel cages**

**2 x 30"H x 30"W stainless steel cages**

Hardwearing and easy to clean the Burtons 'Lifetime' range is the perfect option for everyday caging veterinary practices.

Burtons 'LIFETIME' stainless steel cages have been designed from the ground up to be strong and durable, easily standing up to the rigors of everyday practice life whilst still looking great. Supplied with sound deadening insulation as standard and constructed entirely in stainless steel for easy cleaning and increased hygiene (including the doors, hinges, latches and even screws!)

The cage body is made from the highest grade stainless steel and constructed for superior strength and elimination of dirt traps on

seems. The thicker construction of the Lifetime cages reduces bends and un-shapping, that can cause the thunder noise when the patients move around.

All of our cages are available in standard sizes from stock with special made to measure sizes available on request. The completely modular, standardized design ensures not only an easy set up but allows your bank to grow with your practice. They are even compatible with any existing Burtons, Ken cage & Shorline cages that you may have.

## Features

- Stainless steel, electro-polished removable door may be mounted either side
- Covered rear corners on top and bottom eliminate accumulation of dirt and bacteria
- Stainless steel sub frame makes the front of the cage extra strong
- Radius bends on all sides provides a smooth, safe, easy to clean interior surface
- Stainless steel self-latching door latch is designed for one handed opening and closing, with the ability to accommodate a padlock for greater security
- Stainless steel doorframe and cross members are welded together for strength and rattle-less finish
- Distinctive floor design with fluid retaining lip as standard
- Stainless steel latch keeper latches door securely
- Heavy gauge stainless steel die formed barrel style hinges operate smoothly for the life of the cage
- Latch connector bar is finger notched for smooth, easy opening
- Stainless steel hinge pins remove easily for cage maintenance
- Heavy 12 gauge stainless steel latch brackets are fully adjustable for smooth operation

## Doors

All of our doors utilise heavy stainless steel die formed, adjustable barrel style hinges for years of smooth, reliable operation.

All cages are supplied with electro polished stainless steel mesh doors with vertical spokes spaced at 1" apart and welded at every intersection for superior strength and rattle less operation. The reversible; left or right opening design allows you to easily swap the doors to get the perfect configuration for your practice. All of our doors incorporate radius edges that go corner to corner for easy cleaning and patient security with a single handed opening and a simple push to close self-locking latch design.

All cages are fully compatible with our full range of oxygen therapy, mesh or acrylic doors and can be easily upgraded to one of our specialist doors either at the time of ordering or later as required.

## Socket Layers

Complete your cage arrangements with one of our stainless steel socket layers. Socket layers can add power sockets, optional oxygen supply and temperature control (for heated cages) to your cage bank. All conveniently located at the top of the cage bank.

## Heated Cages

Create a more comfortable temperature controlled environment with the help of built in temperature pads, include safety cut outs as standard and rear mounting on the back of the cage.

Simply add these to any cages at the time of order \*

\*Requires the addition of a socket layer for temperature control panel.

## Bases

All of our cages can be mounted on our purpose made stainless steel mobile, stationary, raised or laminate plinths. Alternatively, you can build your own. Manufacturing details will be sent to you with your order.

### Complete your bank

Complete the look of your bank with a full range of top and side panels. Finished in matching stainless steel, these great looking panels are easy to clean and are the perfect finish to any bank.

### Full range of accessories

We have one of the widest ranges of cage bank essentials from clipboards to IV hooks. Many of which are completed in matching 304 grade stainless steel for a hygienic, easy to clean, long life durability.

### Let us help with the design...

At Burtons, we offer a completely complimentary custom design service, meaning whatever your plan for the room, whether it is to gain ISFM certification or just maximising every bit of space. Our dedicated projects team will use their expert knowledge to complete a design that will be both functional and to specification.

We can arrange for a full site survey or you can simply provide us with your room dimensions and drawings and our projects team will complete a free 3D room design based around your needs.

### Installation

If it's a single cage or a new practice setup the easy installation process is straight forward and fast. Once your happy with your order, all our cages are delivered pre-made and ready to go, but if you prefer why not take advantage of our in-house installation team, who can then arrange a suitable time to build your new bank, working with you to minimize disruption.

(Please note: we recommend that stainless steel cage banks should be 'earth bonded'. You should consult your electrician for advice and have them carry out any recommended work)

### Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	170-028-05
<b>Colour/Style</b>	Stainless Steel

<b>Brand</b>	Burtons
<b>Catagory Sub Heading</b>	Ideal for dog/mixed ward
<b>Type</b>	Bank
<b>Depth</b>	28"
<b>Complete Bank</b>	Yes
<b>Bank Configuration</b>	1x 36"H x 48"W double-door stainless steel cage/ 1x 36"H x 36"W stainless steel cage/ 1x 30"H x 24"W stainless steel cage/ 2x 30"H x 30"W stainless steel cages
<b>Material</b>	Stainless Steel
<b>Mobile</b>	No
<b>Door Type</b>	Full Mesh
<b>Cage Depth</b>	28"

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights