

## PETMAP GRAPHIC 3

129-560-G-3



## PRODUCT DESCRIPTION

---

### **Quick and easy measurement of diastolic blood pressure, mean arterial pressure (MAP) temperature & heart rate, now with a larger 5" (12.7cm) touchscreen display**

The PetMAP graphic range has been specially developed for veterinary use, compact and lightweight the PetMAP graphic offers PPO technology, which allows the measurement of diastolic blood pressure, mean arterial pressure (MAP) temperature & heart rate

The measured values are displayed and trended graphically and numerically on the device. Data trended on the PetMAP can optionally be saved on an SD memory card via the devices integrated SD memory card slot.

The user is able to easily interact with PetMAP Graphic + using its touch screen, with alarms for systolic blood pressure and automatic NIBP monitoring mode.

### **Key Features**

- Display of systolic, diastolic, mean arterial pressure (MAP) along with heart rate
- Correlation to intra-arterial readings
- Easy to use larger 5" touch screen display
- Automatic cuff inflation
- Selectable BP cycle times for measurement reminders
- Graphic display of last 45 readings and NSV
- Tabular display of readings
- Alarm set for systolic blood pressure
- Real-time display of oscillometric
- Correlation to intra-arterial readings
- PPO (petMAP proprietary optimization) - improves the correlation to Intra Arterial Pressures (IAP), the patient's true blood pressure.
- Calculates & automatically displays the normal NSV (nominal session value) after each reading - NSV is most representative

- on the patient BP
- BP cycle time set for automatic reminders to start another BP measurement
- SD memory card to record all BP/HR values and power off Image (SD Card Sold Separately)
- Mini USB Slot
- Supplied with 2 cuff extensions
- Optional pole mount, cover and AC/DC adapter

**Blood Pressure Range & specifications**

- Systolic BP: 40 - 265 mmHg
- MAP: 30 - 200 mmHg
- Diastolic BP: 25 - 195 mmHg
- Heart Rate Range: 40 - 240 bpm
- Temperature: 29C - 46C
- Display: 5" (12.7cm)
- Dimensions: 15 x 8 x 3 cm
- Weight: 460 g
- Batteries: 4 AA NimH rechargeable batteries
- Included cuffs: 2 - 2.5 - 3.5 - 4 - 4.5 - 5.5cm

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	129-560-G-3
<b>Colour/Style</b>	Black
<b>Brand</b>	PetMAP
<b>Monitor Parameters</b>	Heart Rate, Multi-Parameter, NIBP, SpO <sub>2</sub> , Temperature
<b>Screen Type</b>	Touchscreen
<b>Power Source</b>	Battery, Rechargeable
<b>Optional Modules</b>	No
<b>SD Card</b>	Yes

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights