

## WALL MOUNTED O<sup>2</sup> SOLITAIRE ANAESTHESIA MACHINE

127-100-W\*BUNDLE



### PRODUCT DESCRIPTION

## Wall mounted anaesthetic machine from the Burton's Value Range

A low cost alternative to other anaesthetic machines, the solitaire is made from a durable powder coated mild steel with high quality built in 0-10 l/min rotameter and emergency O<sub>2</sub> flush.

The wall-mounted option makes the solitaire extremely convenient in situations where floor space is at a minimum and is an ideal choice for piped gas systems.

#### Key Features

- Wall-mounted powder coated frame
- 0-10 l/min O<sub>2</sub> rotameter with emergency O<sub>2</sub> Flush
- Corner unit outlet (BS22mm Male)
- Product Dimensions 420mm (W) x 460mm (H) x 220mm (D) (see diagram of dimensions above)

**Don't forget your Vaporiser**



The Series 3 is a temperature compensating, out of circuit vaporiser in a traditional design, familiar to most veterinary practitioners. Calibrated specifically for the specified anaesthetic drug using a laser refractometer.

Fitting Type

Choose from the standard taper mounted fitting (cage mount) (UNI) or the quick release (QR) option; allowing you to easily swap the vaporiser to different anaesthetic machines throughout the practice. Additional quick-release compatibility blocks will be required for each anaesthesia machine. (Available separately)

Key filled or Funnel (Pour fill)

Key filled (KF) systems require the addition of key filler adaptors but are generally safer and have a lot less wastage compared to standard funnel filled (FF) vaporisers, saving you money in the long run.

Manufactured to last

We don’t just supply it! It was designed and built by Burtons in the UK at our principal manufacturing facility in Maidstone, Kent Using the latest CNC’s and laser manufacturing processes and high-quality stainless steel materials this product has been engineered from the ground up to last!

*Please note: Pipeline Lead & Vaporiser Sold Separately*

ADDITIONAL INFORMATION

Animal	Mixed
SKU	127-100-w*bundle
Overall Product Dimensions (mm)	420mm (W) x 460mm (H) x 220mm (D)
Gas Type	O <sup>2</sup>
Gas Supply	'E' Size Cylinder (s), 'F' Size Cylinder(s), Piped Gas
Flowmeter Type	Graduated Rotameter
Concentrator	Optional
Optional O2 Alarm	No
Optional N2O	No
Emergency O2 Flush	Yes

### Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights