

ADULT REUSABLE OESOPHAGEAL/RECTAL PROBE (YSI-400) FOR VETRONIC LIGHTNING MONITOR & IMPACT III MONITOR

VSL-W0001B



PRODUCT DESCRIPTION

Adult Reusable Oesophageal/Rectal Probe (YSI-400) for Vetronic Lightning Monitor & Impact III Monitor. Change when broken. Compatible with Lightning & Impact monitors.

ADDITIONAL INFORMATION

Animal	Mixed
SKU	VSL-W0001B
Brand	Vetronic
Accessory Parameter	Temp
Monitor Compatibility	Vetronic Impact, Vetronic Lightning
Change Frequency	When Broken

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

