

ECG LIMB PLATE ELECTRODE, VELCRO STRAP AND SOCKET 4MM BANANA PLUG, PACK OF 4 FOR VETRONIC LIGHTNING MONITOR

VSL-VS0050



PRODUCT DESCRIPTION

ECG Limb Plate Electrode, Velcro strap and socket 4mm banana plug, pack of 4 for Vetronic Lightning Monitor. Change when broken. Compatible with Vetronic Lightning monitors.

ADDITIONAL INFORMATION

Animal	N/A
SKU	VSL-VS0050
Brand	Vetronic
Accessory Parameter	ECG
Monitor Compatibility	BLT E30VET, Edan X12, Mindray ePM12, Mindray uMEC, SurgiVet Advisor
Change Frequency	As Used

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

