

GAS ANALYSER SAMPLING KIT FOR VETRONIC LIGHTNING MONITOR & IMPACT III MONITOR

VSL-VS0057



PRODUCT DESCRIPTION

The VS0057 Gas Analyser Sampling Kit is comprised of:

- **3827** - NomoLine Sampling Line, Adult/Pediatric (5kg and above) **x1**
- **3829** - NomoLine Sampling Line, Infant (1kg-5kg) **x1**
- **VS0055** - x5 RSP 15mm ISO ET Connectors with Sideport(2.5,3,3.5,4,5mm) **x1 set**
- **1948** - Straight Connector 15M - 6mm swivel stem - 15F **x1**

1. For monitors respivet Apnoea monitor

ADDITIONAL INFORMATION

Animal	N/A
SKU	VSL-VS0057
Brand	Vetronic
Accessory Parameter	CO ²
Monitor Compatibility	Vetronic Impact, Vetronic Lightning
Change Frequency	When Broken

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights