

DEVILBISS 10 LITRE OXYGEN CONCENTRATOR

DEVIL-1025UK



PRODUCT DESCRIPTION

The DeVilbiss 10L concentrator has been developed to offer optimal oxygen therapy with a more compact design.

Taking up the same footprint as a 5 litre concentrator and weighing just 19kg compared to other 10L concentrators that weigh around 30kg.

Capable of delivering 87% to 96% of O₂ purity from 2 to 10LPM, the 1025 oxygen concentrator is an ideal investment for your patients and business alike.

It utilises the latest advancements in oxygen therapy and is compatible to be used with cylinder trans fill systems.

Key Features:

- Each system comes with the DeVilbiss OSD (Oxygen sensing device) which monitors the oxygen and ensures patient safety.
- Oxygen outlet incorporating a fire protection adapter
- Visual and audible alarms built into the device indicate - power failure, high gas temperature, high pressure, low/high flow, low oxygen and service required
- Auxillary oxygen port - at the back of the device which can be used to fill cylinder transfill systems
- Small, lightweight and compact design, weighing just 19kg and within the same housing as a 5 litre concentrator
- Durable and reliable device - this unit has been tested under a number of conditions like high temperatures & high humidity
- The unit has a bright LED panel for patient controls, with a protected cannula fitting and recessed humidifier nook to prevent damage
- Lockable flow meter

- Unit accommodates up to 15m of oxygen tubing plus 7ft of cannula

| Specifications | |
|-------------------------|------------------------------------|
| Outlet Pressure | 20.0 ± 1.0 psi (138 kPa +/- 7 kPa) |
| Oxygen Percentage | 2-10 LPM = 87 %-96% |
| Electrical Requirements | 230V, 50 Hz |
| Power Consumption | 670 watts average |
| Weight | 19 kg |
| Dimensions: | 62.2 x 34.2 x 30.4 cm |
| Warranty | 3 years |
| Classification | II Type B |

ADDITIONAL INFORMATION

| | |
|---------------------|--------------|
| Animal | Mixed |
| SKU | DEVIL-1025UK |
| Colour/Style | Black |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights