

## HOTDOG VETERINARY WARMING SYSTEM

119-090-BUNDLE\*



### PRODUCT DESCRIPTION

---

## Keep your patients normothermic with the latest warming technology

**The Hot Dog Veterinary Warming system is unlike no other, using electrically conductive fabric for even heat distribution.**

It reduces hot spots, eradicates air blowing around the surgical site and doesn't cause any noise pollution. It is also one of the most durable, cost effective warming systems available. There are a variety of blanket sizes available, ranging from extra-small to extra-large, so you will be able to fold the blanket over the patient correctly by using the right size.

The touchscreen user interface has been specifically designed to make controlling the patient's temperature a simple procedure. Other features for the user interface include troubleshooting information displayed and alarm error codes to diagnose any alarm readings. The touchscreen user interface allows you to easily access user guides and FAQs, so even in a difficult situation, there is no need to panic, as all the help you need is easily accessible.

You will not hear the controller or blanket when operating, as the system has been designed to run with no loud noises, so you can focus on the important things. It is very lightweight, therefore making it portable, so you can move the system to a safe place when not in use.

#### Key features:

- Touchscreen user interface for easy operation
- Integrated training slide shows
- Troubleshooting information displayed
- Positioning Guidance provided
- Answers to FAQs
- Adjustable settings
- Includes IV-Pole Mount with VESA standard interface for additional mounting options.
- A clear cable retention loop/strap can be used to assist with cable management.
- Quiet operation

- Lightweight
- 1 year warranty (Extended warranty available)
- Controller is compatible with existing cables and blankets.
- Available in multi-port version or single

The average life expectancy of the blankets is 18-24 months. The blankets are also marked with a two-year expiry date label. <http://vetwarming.com/resources> will always have the most current user information.

Watch our introduction video on the HotDog Veterinary Warming System [here](#)

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	119-090-bundle*

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights