


## JELCO I.V CATHETER

400-04-CATHETER



### PRODUCT DESCRIPTION

The Jelco® peripheral IV catheter features FEP polymer construction, “J” point needle design, flashback visualization and an easy to handle length



#### Sharp, beveled needle tip delivers consistent, easy insertions

- Thinwall catheter design provides high flow rates allowing for the use of a smaller gauge device
- FEP polymer catheter and V-point needle design allows reduction of insertion force which minimizes patient discomfort
- Consistent trim accuracy helps to improve vein access

#### Trusted by healthcare professionals for over 30 years

- Jelco® IV Catheters are used in pediatric/neonates for the high flow rates
- ORs use Jelco® IV Catheters for their quick flash, high flow rates, and ease of use for arterial lines and for I.V. access
- A high-performance product from the brand you trust, Jelco®

### ADDITIONAL INFORMATION

Animal

N/A

SKU

400-04-catheter

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights