

QS5 - ULTRAWAVE CLEANING BATH - 4.5L

116-617



PRODUCT DESCRIPTION

Ultrawave's Q-series baths allow accurate cycle parameter control and effective cleaning performance for all applications

Includes a FREE Pack of 6 Ultraclean M2 Pack Detergent

The QS5 is a 4.5Litre heated, digital ultrasonic bath. The menu driven control panel allows simple adjustment of cycle time (0 to 99 minutes), temperature (ambient to 70°C), ultrasonic power and frequency mode for precise and accurate processing.

Combining features such as degas mode, fluid level sensor and a memory card port, the QS5 offers repeatable and traceable cleaning performance with consistent cleaning results every time.

Please note: 3 year warranty included when you order with us*

**1 year as standard, 2 years extra warranty applicable if registered online with Ultrawave within 30 days of placing your order*

| | |
|--------------------------|--|
| Overall Dimensions | 345 x 200 x 265mm (WxDxH) |
| Internal Tank Dimensions | 300 x 150 x 150mm (WxDxH) |
| Basket Dimensions | 275 x 125 x 90mm (WxDxH) |
| Capacity | 4.5Litres |
| Constrution | 304 Stainless Steel Tank |
| Control | Digital Control Time (0 - 99minutes) Temperature (Ambient to 70°C) |
| Operating Features | Memory Card Port For Digital Traceability Water Level Sensor |

| | |
|---------------------|--------------------|
| Ultrasonic Power | 100W |
| Operating Frequency | 32 to 38kHz |
| Heating Power1 | 50W |
| Required Services | 13amp 3 Pin Socket |

ADDITIONAL INFORMATION

| | |
|---------------------------------|-----------|
| Animal | Mixed |
| SKU | 116-617 |
| Brand | Ultrawave |
| Ultrasonic Bath Controls | Digital |
| Ultrasonic Bath Litres | 12.5 |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights