

GIMA DIATERMO MB ELECTROCAUTERY UNIT & ACCESSORIES

GIMA-MB120D-BUNDLE*



PRODUCT DESCRIPTION

The GIMA MB 120 is a high-frequency electrosurgical unit designed specifically for use in veterinary practices

Electrosurgery unit bundle includes Autoclavable MB Handle, Bipolar Cable and 15cm Straight forceps.

Suitable for a wide range of small and medium animal procedures, it supports both monopolar and bipolar techniques, making it a versatile choice for veterinary surgery.

This unit allows for pure cutting (CUT), blended cut/coagulation (BLEND), superficial coagulation (FORCED COAGULATION), and gentle, non-necrotising coagulation (SOFT COAG). Bipolar coagulation is also supported for procedures requiring precision and minimal thermal spread, ideal for delicate soft tissue work in animals.

With robust safety features including constant monitoring of the patient/plate connection, and digital readouts of power output, the GIMA MB 120 ensures reliable and consistent performance.

Key Features:

- Optimised for veterinary use in small animal and exotics surgery
- Minimally invasive procedures with controlled thermal spread
- Compatible with bipolar forceps, ideal for precision veterinary work
- Power and function controls are accessible via the handpiece to maintain focus during surgery
- Last-used settings are stored for procedural efficiency
- Activation options: foot pedal or hand control
- Clear digital display for precise power adjustments
- Adjustment of the acoustic indicator intensity.

Technical Specifications:

- Minimum preselect able power: 0
- Level step: 1
- Maximum output power CUT (W): 120 250 Ohm
- Maximum output power BLEND (W): 90 200 Ohm
- Maximum output power COAG FORCED (W): 80 150 Ohm

- Maximum output power COAG SOFT (W): 60 100 Ohm
- Maximum output power BIPOLAR (W): 40 100 Ohm
- Working frequency: 600 kHz
- Weight kg: 7
- Size LxHxD mm: 260 x 110 x 265
- Operating voltage (switchable): 115 -230
- Mains frequency (Hz): 50-60
- Fuses 230Vac (5x20) TIMED: 2 x 3,15A
- Fuses 115Vac (5x20) TIMED: 2 x 6,3A
- Electrical input power (VA): 300
- Electrical input current (230 Vac) (A): 1,3
- Electrical input current (115 Vac) (A): 2,6
- Five steps adjustable sound level (from 55- to 75dBA): YES
- Self-check: YES
- Power accuracy output warning: YES
- Split or not split patient plate allowed: YES
- Last working condition storing: YES
- Electrical Class (EN60601 - 1): I CF
- MDD 93/42/EEC Class: II B
- Patient circuit: F

ADDITIONAL INFORMATION

Animal	Mixed
SKU	GIMA-MB120D-bundle*
Cut Modes	pure cut, coagulated cut, forced coagulation, soft coagulation, bipolar coagulation

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights