

## MINDRAY DRYLINE WATER II TRAP, ADULT / PAEDIATRIC

IMEC-115-043024-00



### PRODUCT DESCRIPTION

**DRYLINE™ II is a complete water trap subsystem and the ultimate solution to the problems of moisture, humidity and bacteria contamination usually associated with airway gas monitoring.**

The patented design of the DRYLINE™ II water trap ensures virtually unaffected response time in any clinical application for patients over 3kg during anesthesia and intensive care.

The sample volume for the adult dryline water trap is 120ml.

#### Key Features:

- . DRYLINE™ II interfaces directly to AION™ Platinum Multigas Analyzers and can be easily integrated to any design of gas monitors of any make.
- . Protects the gas analyzer from humidity, mucus and bacteria through 0.45 µm hydrophobic bacteria filter
- . Fail-safe protection of the gas analyzer system by additional self-sealing filters
- . Flow optimized gas sample path ensures virtually unaffected response time
- . Semi-disposable water trap for use up to one month
- . Available in versions for adult/pediatric (120 - 200 ml/min) and neonatal applications (70 - 120 ml/min)
- . Fits in receptacle with switches for automatic water trap type identification and pump shut down when water trap is removed
- . Flush or panel mounting options for receptacle
- . Patented design

Mindray Dryline Water II Trap. Change every 2-4 weeks. Compatible with Mindray uMEC & Mindray ePM12 monitors.

### ADDITIONAL INFORMATION

Animal

N/A

<b>SKU</b>	IMEC-115-043024-00
<b>Brand</b>	Mindray
<b>Accessory Parameter</b>	CO <sup>2</sup>
<b>Change Frequency</b>	2-4 weeks

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights