

## PVC E.T TUBES

ET101-2-CONFIG\*



## PRODUCT DESCRIPTION

**Non-toxic medical grade ET tubes made from PVC. Supplied with 15mm connector and Murphey eye.**

### Features:

- Transparent
- X-ray line
- Supplied with 15mm connector.
- Murphey eye
- Available cuffed & un-cuffed
- Cuffed sizes 3.0 to 10.0mm
- Un-cuffed 2.0 to 4.0mm

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	ET101-2-config*
<b>Catagory Sub Heading</b>	Select your type and size
<b>Tube Material</b>	Silicon Coated PVC

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights