

### **EQUINE BLACK RUBBER REBREATHING BAGS**

104-230-BUNDLE



## **PRODUCT DESCRIPTION**

# Ideal for manual patient ventilation or in connection with anaesthesia apparatus

- Equine Specific
- 50mm Thick neck
- Manufactured from black antistatic natural latex rubber
- Complete with tail and ring (excludes open ended version)
- Replace when surface has because sticky, ripped or brittle

### ADDITIONAL INFORMATION

Animal	Equine
SKU	104-230-bundle
Catagory Sub Heading	Range of Sizes Available
MRI Compatible	No

https://burtonsveterinary.com/ Page 1 Sales: 01622 834300 | Service: 01622 834350



#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights