

BURTONS HD VIDEOSCOPE 13MMOD X 3500WL

END-KW1335



PRODUCT DESCRIPTION

Burtons HD Videoscope 13mmOD x 3500mmWL - ideal for mixed practice

Ideal for small animal upper and lower Gi tract and some URT in medium to large breed dogs, aquariums, zoos, URT in equine, guttural pouches, and cystocopy

Key Features:

- 13mm outer diameter
- 3500mm working length
- 2.5 mm biopsy channel
- 4-way angulations 210 degrees up / 100 degrees down
- 90 degree left/right

Supplied in transportation case with:

- Biopsy forceps x 2
- Cleaning brush x 2
- Flushing adapters
- Water bottle
- Spare piston kit
- Spare biopsy cap
- Waterproof cap 3 x 2 viewing and recording screen
- Video cables SD card
- User Manual

ADDITIONAL INFORMATION

SKU	END-KW1335
Brand	Burtons
Animal	Mixed

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights