

**OXYGEN THERAPY SUPPLY KIT**

109-010

**PRODUCT DESCRIPTION****Oxygen Therapy Supply Kit For Piped Oxygen Outlet**

Supplied with 0-15 LPM O<sup>2</sup> flowmeter & 3 meter black feed tubing.

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	109-010
<b>MRI Compatible</b>	No

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

