

**SUNTECH VET30/30E CONTINUOUS BP MONITOR WITH PULSE OX & TEMPERATURE**

SUN-VET30E-CONFIG\*



BURTONS SunTech Range	SUNTECH VET30 Spot Check BP Monitor	SUNTECH VET30 Spot Check & Continuous BP Monitor	SUNTECH VET30E Spot Check BP Monitor + Pulse Ox + Temperature
MOTION TOLERANT MEASUREMENTS	⊙	⊙	⊙
SELECTABLE BP AVERAGING	⊙	⊙	⊙
EDITABLE CLINICAL ALARMS	⊙	⊙	⊙
SMALL ANIMAL & EQUINE MODELS	⊙	⊙	⊙
AVAILABLE IN 3 X COLOURS	⊙	⊙	⊙
BLUETOOTH FOR DATA TRANSFER	⊙	⊙	⊙
PC APPLICATION INCLUDED	⊙	⊙	⊙
BATTERY TYPE	AA ALKALINE	LITHIUM ION	LITHIUM ION
BP SPOT CHECKING	⊙	⊙	⊙
CONTINUOUS BP MONITORING	⊙	⊙	⊙
TEMPERATURE MEASUREMENT	⊙	⊙	⊙
SPO2 MEASUREMENT	⊙	⊙	⊙
Find out more at <a href="http://burtonsveterinary.com">burtonsveterinary.com</a>	WAS £1392 NOW £999 £1099	JUST £1499	JUST £1799



**PRODUCT DESCRIPTION**

**Suntech blood pressure monitor for continuous blood pressure, with pulse oximetry & temperature**

VET30 - For small animal

VET30E - For small animal & equine

Quickly and easily provides reliable vital signs measurements - BP, heart rate, oxygen saturation and temperature.

Combining SunTech Medical’s trusted BP technology with pulse oximetry and temperature monitoring, the SunTech Vet30 Veterinary Monitor gives veterinarians an accurate, portable, and easy to use monitor to perform procedural and continuous non-invasive monitoring on dogs and cats before, during, and after procedures.

Designed for use in emergency assessment, dental procedures, pre-op, and recovery scenarios, the Vet30 quickly and easily provides reliable vital signs measurements. With the Vet30, intervals can be set to automatically monitor and measure systolic and diastolic pressures, heart rate, oxygen saturation and temperature. It is light, portable, and small enough to carry between rooms.

The Vet30 performs motion-tolerant BP measurements with SunTech Medical’s animal-specific Advantage VET BP technology, delivering accurate results with no shaving or gel. Temperature measurements are powered by AccuVet.

The intuitive, touchscreen interface provides users all the necessary features for reviewing, analysing, and transferring vital signs data. Using Bluetooth, data can be transferred from the Vet30 to the trending and reporting PC application that comes standard with the monitor.

Watch our new product spotlight video of the Suntech VET BP Monitor [here](#)

## Key Features

- Three-in-one - Trusted SunTech veterinary BP technology now with pulse oximetry and temperature
- Simple to use - Intuitive, colour, touch screen interface.
- Accurate - BP algorithm was specifically developed for companion animals by experts at a leading vet school.
- Continuous monitoring - Monitoring during procedures is simple with user-defined intervals for BP measurement.
- Portable - Battery-operated device easily moved to different monitoring locations.
- Portable - Battery-powered device can be easily moved between exam rooms.
- Rechargeable - Battery recharged when device is plugged into AC power, reducing disposable battery waste.
- Memory - Up to 960 data points are retained, even when device is powered off. Users can select specific measurements from device memory and perform automated BP averaging
- Data Transfer - Data can be transferred to a data trending and reporting software application for a PC using Bluetooth. Use the PC application to export reports in either PDF or CSV format.

## VET30 (for small animal) package includes:

- Vet30 system monitor
- Protective armour in Flamingo Pink, Peacock Blue, or Tree Frog Green
- Veterinary blood pressure cuffs sizes 1-6
- Patient hose, 4ft (1.2m)
- AccuVet Temperature Probe
- SpO2 Y sensor with clips
- AC Adapter

## VET30E (for small animal & equine) package includes:

- Vet30E system monitor
- Protective armour in Flamingo Pink, Peacock Blue, or Tree Frog Green
- Veterinary blood pressure cuffs sizes 1-8
- Patient hose, 4ft (1.2m)
- AccuVet Temperature Probe
- SpO2 Y sensor with clips
- Reflectance/transflectance SpO2

## Specifications:

- Systolic: 40 - 265 mmHg
- Diastolic: 20 - 200 mmHg
- Mean Arterial Pressure: 27 - 222 mmHg.
- Heart Rate: 25 - 300 BPM
- Device Weight: 2.65lbs (with battery) 1202 grams (with battery)
- Dimensions: 6.25" x 5" x 5.25", 15.9cm x 12.7cm x 13.3cm
- Warranty: 2 years on monitor, 1 year on AccuVet accessories, 90 days on cuffs/hoses

## Temperature & Pulse Ox Specifications:

- Temp range: 0C to 50C
- Temp accuracy: plus, the temperature sensor tolerance
- Temp sensor: YSI 400 compatible
- SpO2 Range: 0-99% functional SpO2
- SpO2 Accuracy:  $\pm 2$  at 70-99% less than 70% is unspecified.
- Pulse rate range: 18-400 BPM
- Pulse rate accuracy:  $\pm 2\%$  or 2 BPM, whichever is greater.

## ADDITIONAL INFORMATION

SKU

SUN-VET30E-config\*

**Category Sub Heading**

Select your colour & animal type

**Monitor Parameters**

IBP, NIBP

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights