

SUNTECH VET30/30E CONTINUOUS BP MONITOR WITH PULSE OX & TEMPERATURE

SUN-VET30E-CONFIG*



PRODUCT DESCRIPTION

Suntech® VET30 Blood Pressure Monitor for continuous blood pressure monitoring, with pulse oximetry & temperature monitoring

VET30 - For Small Animal

VET30E - For Small Animal & Equine

Quickly and easily provides reliable vital signs measurements - Blood pressure, heart rate, oxygen saturation, and temperature. It is especially well suited to small animals, whilst accommodating large animals and the 'E' model version accommodating Equine.

Combining SunTech® Medical's trusted BP technology with pulse oximetry and temperature monitoring, the SunTech® Vet30 Veterinary Monitor gives veterinarians an accurate, portable, and easy-to-use monitor to perform procedural and continuous non-invasive monitoring on dogs and cats before, during, and after procedures. Validated to meet American standards, the Vet30 ensures consistent and repeatable results, even for smaller patients such as cats and small dogs. From the smallest feline to larger canine patients, this monitor excels across the board. Extensive research and white papers back its reliability, making it the preferred choice for accurate blood pressure monitoring in veterinary practices. This monitor sets itself apart from alternative monitors that may struggle to provide consistent readings, particularly in smaller animals.

Designed for use in emergency assessment, dental procedures, pre-op, and recovery scenarios, the Vet30 quickly and easily provides reliable vital signs measurements. With the Vet30, intervals can be set to automatically monitor and measure systolic and diastolic pressures, heart rate, oxygen saturation, and temperature. It is light, portable, and small enough to carry between rooms. Unlike alternative monitors, the Vet30 has demonstrated superior accuracy in studies, with results consistently within 10mmHg of the gold standard for blood pressure measurement in both small and large patients.

Utilising SunTech® Medical's trusted, animal-specific Advantage VET BP technology, the Vet30 performs motion-tolerant BP

measurements, delivering accurate oscillemetric BP measurements and results with no shaving or gel. Temperature measurements are powered by AccuVet. Unlike some alternative devices that may fail to provide a reading in challenging cases, the Vet30 has proven to reliably capture measurements across a wide range of patient sizes and conditions.

The intuitive, touchscreen interface provides users with all the necessary features for reviewing, analysing, and transferring vital signs data. Using Bluetooth, data can be transferred from the Vet30 to the trending and reporting PC application that comes standard with the monitor.

Key Features:

- Three-in-one - Trusted SunTech® veterinary BP technology now with pulse oximetry and temperature.
- Validated & Reliable - Meets American standards, with repeatable, accurate results supported by clinical research. SunTech® motion-tolerant technology has a high success rate for BP readings.
- Great for Small Animals - Proven to deliver consistent results, even in smaller patients such as cats and small dogs.
- Accurate - The BP algorithm was specifically developed for companion animals by experts at a leading vet school.
- Simple to Use - Intuitive, colour, touch screen interface has all necessary features for reviewing and analysing BP data.
- Continuous Monitoring - Monitoring during procedures is simple with user-defined intervals for BP measurement.
- Portable - Battery-operated device easily moved to different monitoring locations.
- Silent - Quietly measures BP, and alarms may be silenced to avoid frightening your patients.
- Fast - Take multiple measurements in a matter of minutes!
- Rechargeable - Battery recharged when device is plugged into AC power, reducing disposable battery waste.
- Memory - Up to 960 data points are retained, even when the device is powered off. Users can select specific measurements from device memory and perform automated BP averaging.
- Data Transfer - Data can be transferred to a data trending and reporting software application for a PC using Bluetooth. Use the PC application to export reports in either PDF or CSV format.

VET30 (for small animal) package includes:

- Vet30 system monitor
- Protective armour in Flamingo Pink, Peacock Blue, or Tree Frog Green
- Veterinary blood pressure cuffs sizes 1-6
- Patient hose, 4ft (1.2m)
- AccuVet Temperature Probe
- SpO2 Y sensor with clips
- AC Adapter

VET30E (for small animal & equine) package includes:

- Vet30E system monitor
- Protective armour in Flamingo Pink, Peacock Blue, or Tree Frog Green
- Veterinary blood pressure cuffs sizes 1-8
- Patient hose, 4ft (1.2m)
- AccuVet Temperature Probe
- SpO2 Y sensor with clips
- Reflectance/transflectance SpO2 sensor
- AC Adapter

ADDITIONAL INFORMATION

SKU	SUN-VET30E-config*
Catagory Sub Heading	Select your colour & animal type

Monitor Parameters

Screen Type

Power Source

Optional Modules

SD Card

NIBP

Full Colour, Touchscreen

Mains, Rechargeable

No

No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights