

## BURTONS VIDEOSCOPIES

END-KW-BUNDLE\*



### PRODUCT DESCRIPTION

## Range of high quality Burton's endoscopes

For a complete setup that allows you to export your findings to a video monitor or laptop, choose from 2 x Burton's HD video processors, and 5 of our Burton's videoscopes.

Available as a full HD DICOM solution, when purchased with VETIT Pro Software.

- END-KW3560** - Burton's HD Videoscope 3.5mmOD x 600mmWL - 1.2mm channel
- END-KW6010** - Burton's HD Videoscope 6.0mmOD x 1000mmWL - 2.0mm channel
- END-KW8515** - Burton's HD Videoscope 8.5mmOD x 1500mmWL - 2.0mm channel
- END-KW9215** - Burton's HD Videoscope 9.2mmOD x 1500mmWL - 2.8mm channel

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	END-KW-bundle*
<b>Brand</b>	Burtons
<b>Biopsy Channel</b>	1.2mm
<b>Working Length (mm)</b>	600

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights