

3ML TERUMO SYRINGE

400-008

**PRODUCT DESCRIPTION****A single use disposable syringe for precise administration of drugs**

- Clear scale markings to enable precise dosage
- Transparent barrel for visualisation
- Finger grips support easy handling
- Firm and definite end-stop
- Siliconisation for optimal plunger movement
- Shelf life of 5 years

***When purchased, this product comes in boxes of 100.**

ADDITIONAL INFORMATION**Animal**

Mixed

SKU

400-008

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights