

BURTONS PROFESSIONAL LARGE PLATFORM VETERINARY SCALES - CLEARANCE

125-920-CL



PRODUCT DESCRIPTION

Robust, high quality, scales, the perfect addition to any practice!

The Burtons professional veterinary scale has a low profile 50mm platform, 100g accuracy and non - slip rubber mat, these scales are comfortable and easy to use.

Superbly designed with its super-size stainless steel base, this scale is ideal to weigh anything from cats to large dogs up to 150kg. The scale features zero, Hold and Tare functions. Animals can be weighed in either kg/g or lb/oz.

It is also easy to store and the indicator comes with a wall bracket keeping it safe and off the floor. The scale can be powered by 6 x AA batteries or mains power, making it mobile if needed.

Key features:

- 1050mm x 505mm Scale Area
- Capacity: 150 KG
- Stainless steel platform
- Easy clean, hygienic Top
- Long control panel lead
- Adjustable feet
- Black non-slip mat included
- Wall mount for control panel included
- Graduations: 100g
- Tare/Hold Functions
- Powered by batteries for use anywhere in the surgery
- AC adaptor supplied for power when batteries are low
- Tare Function - Any unwanted weight can be easily deducted from overall reading

- Hold function - displays a stable weight, even when the animal on the scale is moving
- Robust and built to last

Specifications

- 150kg Capacity
- Mains and Battery
- Uses 6 x AA Batteries (Sold Separately)
- Platform Dimensions: 1093mm x 510mm x 50mm
- Kg/lb Switchable

Clearance products may be non-pristine, for example:

- On occasion some items may show signs of minimal use, which may include minor cosmetic marks which do not affect the functionality in any way.
- Although most of our clearance products are sold in their original manufacturers packaging, there may be occasions, due to packaging damage, where plain white or brown packaging may be used

Clearance units sold with 6 month standard guarantee, unless otherwise specified. This does not affect your statutory rights.

Please see our terms and conditions for more information.

***Image is a representation of product only and does not show actual condition of item**

S/N: 22403983

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	125-920-cl
Overall Product Dimensions (mm)	1093mm x 510mm x 50mm
Colour/Style	Stainless Steel
Brand	Burtons
Accuracy	100g
Capacity	150kg
Ideal For	Canine, Feline

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights