

MACINTOSH LARYNGOSCOPE SET

117-620



PRODUCT DESCRIPTION

Conventional set of Macintosh blades and handle

The Macintosh Laryngoscope Blade is the predominate model of the curved blade types. During laryngoscopy, the tip or beak of the blade is compressed into the angle formed by the base of the tongue and the epiglottis, indirectly raising the epiglottis. Some sizes are available flangeless

Key Features:

- Stainless steel construction and a satin non-glare finish
- Handles are made of chrome plated brass with a knurled grip
- Blades and handles compatible with all major brands of standard illumination laryngoscopes
- Size 1, 2, 3, 4 Macintosh blades

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	117-620

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

