


## HUGER 69 SERIES LARGE ANIMAL ENDOSCOPES

END-HGAGVE-LA-CONFIG\*



### PRODUCT DESCRIPTION

**The 69 series of Huger videoscopes provide a full HD image when used alongside the HUGER 69 series Video processors.**

Using this complete package, you will get an exceptional view of your procedure, with minimised noise and reduced halation. The full HD system gives you a viewing screen like no other. Capture and record images and videos at the click of a button.

Each endoscope is supplied in a hard shell transportation case and includes factory accessories.

#### Key Features

- Get a full-HD image resolution when used with the Huger 69 processors
- Tip Deflection: Up/Down:180 °/180 °, Left/Right:160 °/160 °
- Field of View: 145 °

### ADDITIONAL INFORMATION

<b>Animal</b>	Equine
<b>SKU</b>	END-HGAGVE-LA-Config*
<b>Working Length (mm)</b>	1050mm
<b>Outer Diameter (mm)</b>	6.0mm, 8.0mm, 9.2mm, 9.6mm, 12.9mm

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights