


HUGER 69 SERIES LARGE ANIMAL ENDOSCOPES

END-HGAGVE-LA-CONFIG*



PRODUCT DESCRIPTION

The 69 series of Huger videoscopes provide a full HD image when used alongside the HUGER 69 series Video processors.

Using this complete package, you will get an exceptional view of your procedure, with minimised noise and reduced halation. The full HD system gives you a viewing screen like no other. Capture and record images and videos at the click of a button.

Each endoscope is supplied in a hard shell transportation case and includes factory accessories.

Key Features

- Get a full-HD image resolution when used with the Huger 69 processors
- Tip Deflection: Up/Down:180 °/180 °, Left/Right:160 °/160 °
- Field of View: 145 °

ADDITIONAL INFORMATION

Animal	Equine
SKU	END-HGAGVE-LA-Config*
Procedure Type	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy
Working Length (mm)	1050mm
Outer Diameter (mm)	6.0mm, 8.0mm, 9.2mm, 9.6mm, 12.9mm

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights