

## BURTONS K9 GASTROSCOPY TRAINING MODEL

END-TRM1



### PRODUCT DESCRIPTION

---

**The Burton's K9 Gastroscopy Training Model has been designed based on a 25-kilo dog, and contains a range of clinical disease states and pathology for veterinarians to investigate.**

Made of synthetic material, the model is designed for flexible veterinary endoscopes up to 9.8mm in outer diameter.

Zippered entry enables the introduction/ removal of foreign bodies for retrieval training.

Ideal for Universities, Colleges and Veterinary Endoscopy Training centres, the model is supplied in a hard case with lubricant and can be used with/without the case for teaching purposes.

#### **Anatomically the model is made up of:**

- the oesophagus
- lower oesophageal sphincter
- fundus
- cardia
- body
- antrum
- pyloric sphincter
- duodenum

#### **Disease states include:**

- ulceration of the lesser curvature and thickened area
- petechial haemorrhage pyloric antrum
- peyers patches
- duodenal papilla
- bile stains
- nodular changes in duodenum

## ADDITIONAL INFORMATION

<b>SKU</b>	END-TRM1
<b>Brand</b>	Burtons
<b>Animal</b>	N/A

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights