

## VETERINARY CLEAR PERSPEX MASKS WITH DIAPHRAGM

104-223-BUNDLE



### PRODUCT DESCRIPTION

## Polycarbonate high visibility mask with soft silicone rubber diaphragm to ensure a good seal

These transparent masks are manufactured from high impact plastic and are fitted with a soft flexible rubber diaphragm, which makes a good seal to the patient and therefore keeps the leak of anaesthetic gases to a minimum, with largest sizes supplied with an adjustable harness. We offer a wide range of sizes to suit most applications ranging from mice, ferrets, cats, dogs to larger animals such as sheep and goats.

#### Key Features:

- Replacement diaphragms available
- Adjustable harnesses available on larger masks
- Available individually or as a set of 6

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	104-223-bundle
<b>Catagory Sub Heading</b>	Range of Sizes Available

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights