


**AG CUFFILL**

104-000

**PRODUCT DESCRIPTION**

**The AG Cuffill is a cuff inflator and manometer, suitable for use with several airway devices; ETTs, tracheostomy tubes, and laryngeal mask airways, including supraglottic airway devices**

It will inflate and monitor the cuff pressure of the airway device between the ranges of 0-99cmH<sub>2</sub>O, with an accuracy of  $\pm 2$ cmH<sub>2</sub>O, and is designed to be used up to 100 times. The remaining amount of uses blinks onto the display before displaying "00", indicating it is ready for use.

Not suitable for checking the pressure of liquid-based balloons, such as those on foley urinary catheters.

**ADDITIONAL INFORMATION**

<b>SKU</b>	104-000
<b>Animal</b>	Mixed

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights