

## HOT DOG MEDIUM 406MM X 584MM BLANKET - CLEARANCE

119-084-CL



## PRODUCT DESCRIPTION

### Additional warming blanket for Hot Dog System

Designed for the HotDog® warming system these blankets use breakthrough technology, which reduce these risks designed to reduce the risk of hypothermia. The patented ThermAssure™ fabric is a conductive polymer which responds to low-voltage current, warming throughout, without hot spots.

Choose from a large range of reusable mats and single or dual port controllers options.

#### Includes:

- Medium blanket (406mm x 584mm)
- Single Port Controller with UK Plug
- Large blanket (558mm x 787mm)

#### Key Features:

- Fast warm up time - warms up in minutes
- Sensor ensures perfect temperature continuously
- Flexible enough to use under, over or around patient
- Heavy duty nylon outer shell is easy to clean and resists damage/punctures
- Reusable - Economical and saves money
- Adjustable straps to ensure maximum contact are included
- Even distribution of heat
- Durable waterproof blankets
- Maintenance free controller
- Quiet Operation
- No blowing of air on surgical site
- Easily fits through cages slats
- Wont incubate bacteria!

**Safety Features:**

A sensor within the blanket ensures the blanket stays at the perfect temperature at all times, a second sensor works as a safety back-up. The controller has 2 over - temperature alarms for safety. If the blanket exceeds maximum temperature setting, it will set off an alarm with a constant beep and light. If the temperature continues to rise, it will automatically shut off ensuring no patient can overheat. If the alarm sounds, the blanket can be reset to safe temperatures by turning the controller off and on.

**S/N: 2040575****ADDITIONAL INFORMATION**

<b>SKU</b>	119-084-cl
<b>Animal</b>	Small Animal
<b>Generic Size</b>	Medium
<b>Colour/Style</b>	Blue

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights