

REFRACTOMETER FOR URINE AND SERUM SPECIFIC GRAVITY

117-630



PRODUCT DESCRIPTION

A compact, handheld diagnostic device for the simple and accurate measurement of serum/plasma protein and urine specific gravity samples.

Key Features:

- Large rubber handle for grip and comfort
- Rubber eyepiece for comfort
- Focusing and wide field of view
- Blue coloured prism with 3 scales: protein/S. G/refractive index and sample cover lid
- Temperature compensation with calibration
- Supplied with storage case

ADDITIONAL INFORMATION

| | |
|---------------------|--------------|
| Animal | Small Animal |
| SKU | 117-630 |
| Generic Size | X-Small |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights