

**'C' MOUNT FOR 1/2IN CHIP SENSOR 0.2A 0.6X CAMERA ADAPTOR FOR MOTIC TRINOCULAR**

145-202

**PRODUCT DESCRIPTION**

Complete your digital microscopy kit with the Motic C-Mount Adapter

**Key Features**

- C-Mount 1/2in Chip Sensor Adaptor - 0.2A 0.6x
- Required to connect Motic Camera to Microscope
- Compatible with Motic SMZ 171, BA210, 310, 410 and AE2000 series microscopes

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	145-202
<b>Microscope Complete Kit</b>	No

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights