

OXI-GEN10 10 LPM OXYGEN CONCENTRATOR - CLEARANCE

112-404-CL-2



PRODUCT DESCRIPTION

Produce the same volume of Oxygen in an 'E' size cylinder for just 10p!

The Oxy-Gen 10 is a completely self-contained mobile oxygen concentrator producing a constant flow of medical grade (93%) oxygen at 10 l/pm.

Designed for 24hr operation with a built in purity meter and ultra-quiet operation Oxy-Gen10 concentrators are being used more and more as an alternative to using high pressure and costly cylinder oxygen. Other benefits include reducing manual handling and fire risk in your practice, as well as ensuring a reliable supply when needed.

- Oxygen Output: Max. up to 10 l/pm @ 5.8 to 11.6psi
- Oxygen Purity: Max.93% +/-3% (dependant on flow rate)
- Sound Levels:<50dba@ 1 metre
- Dimensions:43 x 37 x 65cm WDH
- Weight:30kg
- Power:220v AC 50/60Hz single phase, 480W

S/N: 056543

ADDITIONAL INFORMATION

Animal

Mixed

SKU	112-404-cl-2
Overall Product Dimensions (mm)	430 x 370 x 650 mm WDH
Colour/Style	White
Gas Type	O ²
Concentrator	Yes
Optional O2 Alarm	No
Optional N2O	No
Emergency O2 Flush	No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights