

## HEINIGER STYLE MINI TRIMMER

HNG-707-281



### PRODUCT DESCRIPTION

**Suitable for dogs and cats, the Heiniger Style Mini Trimmer gives you high levels of performance and precision.**

The Style Mini allows you to shear with extreme precision with incredibly low noise levels, making it ideal for blood sampling and cannula placements.

The adjustable snap-on comb ensures a cutting length ranging from 3 to 6 mm. This trimmer is particularly suitable for final styling. After use, the shearing head can be easily cleaned.

This trimmer has a powerful motor, letting you glide through any patient with ease. It is light with a weight of just 120g, and the high performance Li-Ion battery gives you a long lasting battery life of 120 minutes. Vibrations are minimised when using the Style Mini, and when using it, sound levels are always at a minimum.

We recommend lubricating the clipper head of your speed clippers, as insufficient lubrication results in overheating and reduced durability of the blades!

To keep your clipper blades in top condition, simply:

- Oil the clipper blades with a few drops during and after clipping, and wipe off any excess oil with a dry cloth.
- Check which clipper you have, as some run at higher speeds, meaning the blades will heat up quicker. Regular lubrication is even more important for these.

#### Specifications:

- Operating time: up to 120 minutes

- Charging time: 90-120 minutes
- Cutting length: 0.3 mm
- Adjustable comb attachment: 3-6 mm
- Speed: 5500 dbs/min
- Dimensions: 146L x 30W x 38H mm
- Weight: 120 g

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	HNG-707-281
<b>Clipper Speed(s)</b>	5500 dbs/min
<b>Number of Clipper Speeds</b>	1
<b>battery run time (also put charge time)</b>	120 mins
<b>Compatible with other brand blades ?</b>	No
<b>Replaceable Battery</b>	N/A

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights