

## HEINIGER STYLE MINI TRIMMER

HNG-707-281



### PRODUCT DESCRIPTION

## Suitable for dogs and cats, the Heiniger Style Mini Trimmer gives you high levels of performance and precision.

The Style Mini allows you to shear with extreme precision with incredibly low noise levels, making it ideal for blood sampling and cannula placements.

The adjustable snap-on comb ensures a cutting length ranging from 3 to 6 mm. This trimmer is particularly suitable for final styling. After use, the shearing head can be easily cleaned.

This trimmer has a powerful motor, letting you glide through any patient with ease. It is light with a weight of just 120g, and the high performance Li-Ion battery gives you a long lasting battery life of 120 minutes. Vibrations are minimised when using the Style Mini, and when using it, sound levels are always at a minimum.

We recommend lubricating the clipper head of your speed clippers, as insufficient lubrication results in overheating and reduced durability of the blades!

To keep your clipper blades in top condition, simply:

- Oil the clipper blades with a few drops during and after clipping, and wipe off any excess oil with a dry cloth.
- Check which clipper you have, as some run at higher speeds, meaning the blades will heat up quicker. Regular lubrication is even more important for these.

#### Specifications:

- Operating time: up to 120 minutes

- Charging time: 90-120 minutes
- Cutting length: 0.3 mm
- Adjustable comb attachment: 3-6 mm
- Speed: 5500 dbs/min
- Dimensions: 146L x 30W x 38H mm
- Weight: 120 g

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	HNG-707-281
<b>Clipper Speed(s)</b>	5500 dbs/min
<b>Number of Clipper Speeds</b>	1
<b>battery run time (also put charge time)</b>	120 mins
<b>Compatible with other brand blades ?</b>	No
<b>Replaceable Battery</b>	N/A
<b>Blade Size</b>	N/A

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights