

EDAN ROLLING STAND ADJUSTABLE WITH BRAKES

E-83.60.262078



PRODUCT DESCRIPTION

The Edan rolling stand is extremely durable, controllable and stable

It has a sturdy basket for storage of medical equipment/accessories and the stand is height adjustable (between 90cm and 109cm). The brakes allow you to keep it in place, within a crowded and busy environment.

Features

- Stainless metal design
- Anti-scratch protection
- Height adjustable
- Lockable casters

Specifications

- Size: 530mm (L) x 530mm (W) x 900mm - 1100mm (H)
- Weight: 10 Kg
- Caster diameter: 75mm

Edan Rolling Stand Adjustable with brakes. Change when broken. Compatible with Edan monitors.

ADDITIONAL INFORMATION

Animal	N/A
SKU	E-83.60.262078
Accessory Parameter	Other
Monitor Compatibility	Edan X12
Change Frequency	When Broken

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights