

MINDRAY UMEC12 MULTI-PARAMETER MONITOR WITH TROLLEY & EXHAUST KIT

UMEC12VET-CO2-PACKAGE



PRODUCT DESCRIPTION

Safety & Efficiency for Veterinary Monitoring

The Mindray uMEC12 Veterinary Monitor is a small, lightweight, and compact multi-parameter monitor packed with features.

Bundle includes trolley and CO₂ exhaust kit

Parameters measured: CO₂/Temp/ECG/SPO₂/Respiration - with Blood Pressure cuffs.

This neatly packaged multi-parameter monitor has a 12.1inch high-resolution touchscreen, is lightweight, and has up to 4 hours of continuous monitoring and 120 hours trend / 48 hours waveforms reviewing.

This intuitive monitor has been designed with your practice needs in mind. The following features help to ensure the efficiency of monitoring performance. There is no fan built in, so you do not need to worry about any animal hair getting tangled. The waterproof casing will give you peace of mind that your monitor will not be affected if there are some unexpected spillages, and accessory storage at the rear of the monitor keeps all of your included accessories tidy and to hand.

Designed with ease of use in mind, the Mindray uMEC12 patient monitor helps streamline workflows and improve efficiency through its responsive touchscreen display. Its intuitive user interface allows for faster, simpler operation, even for new users, reducing training time and enabling teams to focus more on patient care. Compared to earlier models such as the IMEC8 and IPM12, the uMEC12 features an angled screen, offering improved visibility and more comfortable viewing for the operator.

Watch our product spotlight video on the uMEC12 Multi-Parameter Monitor [here](#)

Key Features

- 12.1-inch high-res LED screen with touchscreen
- Unique accessory cabinet to keep everything organised
- Less than 4kg weight with battery for extra portability
- Low power design
- Durable and robust casing
- Sidestream CO2 - Gold-standard artema technology for accurate measurements
- Bend-resistant accessories

Advanced Technologies:

- Mindray's patented multi-lead ECG Algorithm greatly improves the accuracy of measurement and reduces false alarms
- NIBP quick-measurement technique reduces the discomfort caused by cuff inflation, especially for patients suffering from hypertension or hypotension
- The anti-interference SpO2 algorithm provides accurate measurement even when the patient is mobile
- Large capacity for data storage enables a comprehensive review of patients' historical data, and external USB storage devices are also supported
- 8-hour continuous runtime with one Lithium-ion battery
- 5-lead ECG cable allows for full screen ECG mode with augmented traces. (Additional extra)

Box includes:

- Power cable
- Battery pack
- Dryline airway adapter
- Dryline II Water Trap, adult
- Neonatal NIBP hose
- 3 Lead ECG set
- Temperature probe
- Disposable cuffs
- Reusable SPO2 sensor
- Sampling line
- Water trap
- Airway adapter
- CO2 nasal sample cannulas
- Crocodile clips

ADDITIONAL INFORMATION

Animal	Mixed
SKU	UMEC12VET-CO2-Package
Colour/Style	White
Brand	Mindray
Catagory Sub Heading	Includes Trolley & Exhaust Kit
Screen Type	Touchscreen
Power Source	Mains, Rechargeable
Optional Modules	No
SD Card	No

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights