

## EDAN X12 MULTI-PARAMETER MONITOR WITH TROLLEY & EXHAUST KIT

E-X12VETCO2-BUNDLE



### PRODUCT DESCRIPTION

## The perfect everyday monitoring solution for veterinary practices.

#### Bundle includes:

- Edan X12 Veterinary Multi Parameter Monitor CO2
- Edan Rolling Stand
- CO2 Exhaust Kit

The Edan X12 Multi-Parameter is an ultra-slim, lightweight patient monitor. Designed to meet the demands of modern veterinary practices, this monitor features a crystal-clear, tilted touchscreen display, ensuring optimal visibility from all angles. Its highly responsive touch interface facilitates swift response for quick and easy usability. A comprehensive, yet uncomplicated, dashboard offers comprehensive data visualization options, enabling customizable range settings and alarms, while a simple night mode minimizes disruption during patient monitoring.

Furthermore, the X12 has been developed with a focus on user-friendliness, making it easy to use and understand for all practitioners, reducing training time. This innovative device empowers veterinary professionals to dedicate more time to patient care.

The Edan X12 Multi-Parameter Monitor covers a wide array of vital parameters, including ECG, NIBP, SPO2, Temperature, and Sidestream CO2. Setting it apart from most small animal multiparameter monitors, the Edan X12 offers a low sidestream capnography sampling rate of just 50ml/minute - providing reliable capnography traces and end-tidal carbon dioxide readings for your smallest patients - down to 500g.

Additionally, you won't need to swap the water trap to get different side stream sampling rates as there is only one size with the Edan X12, and you can change the sampling rate in the monitor instead.

This intuitive monitor has been designed with your practice needs in mind. The following features help to optimize monitoring efficiency. There is no built-in fan, to prevent debris accumulation and battery drainage. Its waterproof casing will give you peace of mind that your monitor will not be affected by unexpected spillages, whilst an impressive battery life of up to 8 hours for uninterrupted operation. Additionally, the monitor's clearly colour-coded ports make attaching accessories simple. It's ideal for any veterinary facility looking for a reliable, lightweight patient monitor.

All required accessories are included with the monitor for you to start monitoring as soon as it arrives. A wide range of additional accessories and consumables are readily available. Including an optional trolley that enables easy manoeuvring around your patients and your practice. The trolley is incredibly sturdy, with a 5-spoke metal base, adjustable height (with anti-drop mechanism) and a basket with hooks for spare accessories. A wall mount option is also available.

**Key Features:**

- Clear 12.1-inch, high-resolution, LED, tilted touchscreen
- A battery life of up to 8 hours
- Durable and robust casing
- Ultra slim, compact design.
- Low sampling rate option of just 50ml! Perfect for small furries and cats as well as the standard higher rates for larger patients.
- Fanless design - reducing dirt and hair traps, lowering fault risks, and reducing battery consumption.
- Colour-coded, easily accessible accessory ports
- A large range of high-quality accessories are available, including an oesophageal ECG and temperature probe.
- Perfect for everyday practice monitoring of companion animals.
- Mounting options available, including a height-adjustable trolley and wall mount bracket.

**CO<sub>2</sub> (Sidestream)**

- Sampling rate as low as 50ml/min
- Superior water trap design, preventing monitor from humidity and damage

**ECG + Temp**

- Auto lead placement detection to reduce false alarms
- Customisable filter for accurate ECG waves and measurements
- Auto-detection of 33 types of arrhythmias
- Additional oesophageal probes that monitor temperature + ECG
- 2 x temp probe ports, allowing you to measure core temp and peripheral temp simultaneously

**NIBP**

- Dual dust filter design to reduce blockages
- Cleaning mode for routine maintenance

**SpO<sub>2</sub>**

- Reference reading of Perfusion Index (PI) according to perfusion changes
- Improved shape design of sensor for stable SpO<sub>2</sub> readings
- The Edan X12 also offers a long-lasting battery, 240 hours of trend review, and built-in WiFi.

**Included in the box is a full set of accessories to get you started, including:**

- 2-part 3 x lead ECG cable - press stud style
- Box of 5, high quality, atraumatic crocodile clips that can be used with press studs - great for conscious patients
- SpO<sub>2</sub> sensor with 2 x crocodile clips included
- Water trap
- 1 x sample line
- 1 x nasal line
- 1 x temperature probe
- Battery pack - also includes foam pad to hold the battery pack in position
- Blood Pressure Cuffs
- Screwdriver for fitting battery pack
- Comprehensive user manual
- Shortcut flashcards

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	E-X12VETCO2-bundle
<b>Colour/Style</b>	Black & Grey
<b>Monitor Parameters</b>	CO <sup>2</sup> , ECG, Heart Rate, NIBP, Respiration , SpO <sup>2</sup> , Temperature
<b>Screen Type</b>	Full Colour, Touchscreen
<b>Power Source</b>	Battery

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights