

EDAN AX3 VET ULTRASOUND MACHINE

USEN-AX3 VET-BUNDLE*



PRODUCT DESCRIPTION

Designed to deliver uncompromised performance at a cost-effective price, The Acclarix AX3 VET Ultrasound Machine provides reliability alongside high specifications and additional features.

The high spec resolution paired with TAI (Tissue Adaptive Imaging) results in superb, detailed imagery, particularly at depth, by continuously and automatically optimising the image. With a thoroughly thought-out design, the AX3 VET improves workflow efficiency whilst giving the clearest possible results. The Dual transducer ports will simplify workflow, whilst an extra screen above the keyboard provides a highly sensitive touch screen for easy setting modifications at your fingertips. All whilst still being ultra-light and compact for easy portability. The intelligent software provides dedicated pre-sets and measurement packages for veterinary applications, with a customisable user interface for personalised design, and multiple data management methods: DICOM, network, PDF, AVI, and BMP file export.

This machine is supplied with 2 probes, Linear L12-5Q (5-12MHz) & Micro-convex (5-10MHz)

Key Features:

- Superb, detailed resolution, particularly at depth
- Tissue Adaptive Imaging (TAI) continuously and automatically optimises imaging, allowing more focus on different kinds of animals.
- High-sensitive touch screen for daily operation
- Dedicated pre-sets and measurement packages for veterinary applications
- Customisable user interface for personalised design

- Multiple data management methods: DICOM, network, PDF, AVI, and BMP file export.
- Dual transducer design to improve diagnostic efficiency
- Unique dual-batteries design supports long continuous scanning, more than 2 hours
- Ergonomic control panel with sensitive touchscreen.
- Ultra-light and ultra-compact design for easy portability

ADDITIONAL INFORMATION

Animal	Mixed
SKU	USEN-AX3 VET-BUNDLE*
Colour/Style	Light Grey
Ultrasound screen size	15 inch
Ultrasound colour	Colour

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights