

## EDAN U60 LUXURY MOBILE TROLLEY

USEN-E-83.63.560344



### PRODUCT DESCRIPTION

**The Edan U60 Luxury Mobile Trolley provides a comfortable, multi-functional mounting solution for your ultrasound systems.**

The trolley comes with a printer tray and four wheels with clips for steady placement. It also features a built-in probe holder and an accessories drawer.

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	USEN-E-83.63.560344
<b>Colour/Style</b>	Light Grey

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

