

EDAN AX3 HEIGHT ADJUSTABLE TROLLEY

USEN-E-83.63.560585



PRODUCT DESCRIPTION

The Edan AX3 Height Adjustable Trolley is a simple, ergonomic and cost-effective solution designed for the Acclarix AX3 Ultrasound System.

With its compact design, the Trolley can fit into small spaces and can easily be moved around the practice. It is height adjustable and has a fixed deck angle of 15° for easy viewing of touch screens. It also features 4 transducer holders, 2 gel holders and a tray for comfortable placement of tools.

ADDITIONAL INFORMATION

Animal	N/A
SKU	USEN-E-83.63.560585
Colour/Style	Light Grey

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights