

600 GIVING SETS + FREE BURTONS EN-V5+ INFUSION PUMP

EN-V5-PROMO-600



PRODUCT DESCRIPTION

One of the only infusion pumps on the market that allows storage of multiple giving set calibrations

Bundle Includes:

- Free Burton's EN-V5+ Infusion Pump - FREE!
- 500 x Burton's SuperBlue Standard Giving Sets
- 100 x Burton's Burton's SuperBlue Giving Sets
- IV Line Sample Pack - FREE!

The Burton's EN-V5 sets the new standard for modern, user intuitive infusion pumps. It is one of the only pumps on the market to allow you to store multiple giving set calibrations, meaning you no longer need to buy separate pumps for paediatric, standard infusion and bloods. Store up to 200 calibrations of giving sets.

The EN-V5's large touchscreen really adds to how user intuitive it is, easily change settings and modes as well as clearly view what the pump is currently doing, including a 'Time remaining' box that allows you to judge when the pump needs checking next.

Its battery will last up to 10 hours when fully charged, allowing you to move it easily between theatre and recovery, a quick and easy Allen key fitting allows you to change over from the infusion pole bracket to the cage mount, both of which are included with the pump.

Your Burton's EN-V5 will come already calibrated to Burton's Giving Sets, please let us know when you place your order if you would like it calibrated to a different type.

The exclusive Burton's edition kit comes with a pole mount, cage mount bracket and 'L' shaped power lead as standard - The perfect set up for practices!

Watch our product spotlight video on the EN-V5 Infusion Pump [here](#)

Key Features:

- Stores up to 200 giving set calibrations
- Interchangeable IV pole bracket, cage bracket and 'L' shaped power cord included as standard
- Night mode - dims brightness of the screen so it can still be seen clearly but doesn't shine bright lights into a dark room, causes less stress for recovering patients.
- Suitable for bloods
- Right angle power lead - keeps leads away from prying paws in cages.
- Parameter settings and editing directly on the touchscreen
- 5 infusion modes + micro mode control
- Infusion rates range from 0.1ml/h to 2000ml/h
- IP44 splashproof
- Automatic anti-bolus
- Screen lock to avoid misoperation
- Purge mode for removing air bubbles

Specifications:

- Case material: PBT+PC
- Infusion Modes: Rate mode, Time mode, Body weight, Drip sequence, Drug library & Micro
- Accuracy +/- 5%
- VTBI: 0 - 9999.99 ml
- KVO: 0 -5l/h adjustable (min increment 0.01ml/h)
- Bolus rate: 0.1-2000ml/h (automatic & manual bolus available)
- Occlusion level: 5 levels selectable (50mmHg - 900mmHg)
- Alarm type: VTBI infused, occlusion, pressure drop, battery low/ empty, KVO finished, standby time end, door open, single air bubble, accumulated air bubble, abnormal drip rate, empty drip chamber, VTBI near end, reminder alarm, no power, drop sensor connection, system error
- Weight - 1.5kg

ADDITIONAL INFORMATION

Animal	N/A
SKU	EN-V5-PROMO-600
Battery	Yes
Screen Type	Digital Display, Full Colour, Touchscreen
Volume Control	Yes

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

