

**ANTI-STATIC O<sup>2</sup> FEED TUBING**

104-19\*BUNDLE

**PRODUCT DESCRIPTION****Replacement anti-static O<sup>2</sup> tubing**

Available in:

- Internal - 6.0mm External - 12.0mm
- Internal - 5.5mm External - 10.0mm
- Sold 1 meter lengths (Quantity selected with be supplied uncut)

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	104-19*bundle
<b>Category Sub Heading</b>	Range of Sizes Available

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

