

TAFONIUS JUNIOR LARGE ANIMAL VENTILATOR

VSL-TAF-JNR-BUNDLE*



PRODUCT DESCRIPTION

The Tafonius Junior is not another bag-in-a-bottle ventilator

You won't find an ascending or descending bellows anywhere on this machine. In fact Tafonius Junior behaves more like a normal re-breathing bag than anything else. During spontaneous respiration the cylinder volume rises and falls just like a re-breathing bag. During IPPV the cylinder volume decreases during inspiration exactly like a bag being squeezed and refills in expiration like a re-breathing bag. The way that the 'bag' is squeezed is entirely under your control. Also the emptying and filling of this 'bag' requires no effort from the patient because the motor controlling the cylinder does all the work. This is what we call the Tafonius Airway Servo System.

Benefits

It costs around 85p/hour* to drive a ventilator with oxygen. Over a year of operating at only 6 hours per day that mounts up to over £1300 in oxygen costs alone. *(based on a 600kg horse at 6 breaths per minute and a venturi efficiency of 40%) designed drainage channels that allow you to easily flood the piston area with cleaning agent. The soda-lime absorber is readily removed and the piston is automatically retracted when the cylinder is opened to give full and easy access to the 'inside' of the patient circuit.

Key Features

- Saves Money - Many hospitals use oxygen as their driving gas. A rough guide for a 500kg
- Can be attached to existing ventilator
- Small footprint and lower height for greater manoeuvrability
- Can run off battery - No need for long extension leads
- Easy cleaning - Tafonius Junior has readily accessible airways for cleaning and specially
- 20L tidal volume
- No Pressurised driving gas required
- Suitable for foals
- Silent operation
- Battery backup
- CPAP/PEEP and ASSIST Modes
- Two versions: with and without anaesthetic circuit

- Easy-clean : Cylinder rotates for easy access
- Small footprint
- Uses the same technology as the original Tafonius

ADDITIONAL INFORMATION

Animal	N/A
SKU	VSL-TAF-JNR-BUNDLE*
Colour/Style	White
Brand	Vetronic Services Ltd
Gas Type	Various Options
Gas Supply	'E' Size Cylinder (s), Piped Gas
Flowmeter Type	Graduated Rotameter
Concentrator	No
Optional O2 Alarm	No
Optional N2O	Yes
Emergency O2 Flush	Yes

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights