

COCOON TROLLEY WITH BASKET

CWS109



PRODUCT DESCRIPTION

Trolley with basket for the Cocoon Warming System

- Integrated basket for holding spare blankets and cables
- Ideal for theatre and recovery wards
- Can be used in conjunction with the Burtons Cocoon Cage Bracket
- 670mmH with 5 anti-static castors
- Weight: 2.25kg

ADDITIONAL INFORMATION

SKU	CWS109
Animal	Small Animal

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

