

**MULTIPLE 'E' SIZE CYLINDER TROLLEY**

115-088

**PRODUCT DESCRIPTION**

**Constructed from mild steel with a powder coated finish these trolleys allow the safe and easy moving of up to 12 'E' size cylinders**

- 4 Wheel design
- Individual , 2 level holders
- 'E' Size
- Capacity for up to 12 cylinders

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	115-088
<b>Colour/Style</b>	Blue
<b>Cylinder Size</b>	E

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights