

EXAMION EQUINE DR PORTABLE CASE SYSTEM WITH MEDIUM PANEL

XR-EXSY000007-PROMO24



PRODUCT DESCRIPTION

A completely wireless, two-component DR solution for horses and mixed animal use.

Supplied in one compact, easy-to-transport case, the Examion Portable M Case is perfect for those out in the field visits where you need to capture the images quickly, yet accurately.

With the AED, DR panel, the system can be ready to use within seconds, and can be used with any generator available, the AED allow you to simply position the generator in line with the AED sensor ready for capture. Full images are viewable within 5 seconds, allowing you to ensure you have the images you require before leaving the patient's side. This dramatically reduces the number of re-takes needed, and therefore minimises the stress for the patient and allows for quicker analysis of the images.

The M size of the panel is based on the standard x-ray film format of 24cm x 30cm, making it easily portable and easy to handle, weighing in at just 1.9kg.

This package includes as standard:

- Robust, protective, transportation case
- Soft carrying case
- Medium Wi-Fi DR Panel
- 2 x rechargeable panel batteries
- Compact battery charging station
- Lead back scatter: 0.15mm 24*30
- Cowboy box for Medium Panel - with replaceable handle, telescopic distance holder and telescopic foot
- High-end portable laptop
- AQS software included as standard with the following modules: PACS Vet for image storage, 1 x AQS license, X-ray module, Image optimisation module, METAGRID Module - giving you a virtual module, saving you having to purchase extra physical grids, Multi-organ module - allowing each user to have their own organ capture order pre-sets, depending on how they like to work.

This package also currently includes the following promotional extras:

- AQS Horse Module - horse default settings pre-loaded on turning the system on
- AQS DropBox Cloud for easily storing and sharing patient images.
- AQS Synchronisation module - enabling the panel and laptop to connect directly with one another and automatically sync.
- 5-year panel warranty

Key Features:

- Compact and lightweight, for faster and easier handling around horses
- A small, robust storage case keeps the complete setup safe, whilst easily portable
- High-resolution images due to its small pixel size of just 76 microns.
- The detector and laptop communicate directly meaning this is all you require for a complete examination.
- System is ready to use within seconds of being turned on.
- AED (Automatic Exposure Detection) Panel allows the panel to be used instantly with any generator available.
- AQS software as standard - enables easy x-ray image diagnostics via the onboard radiology viewer, with diagnostic tools and measurement functions.
- Integrated image archiving with direct access to all archived X-ray images and patient data
- Share X-rays with clients and colleagues via the integrated AQS Dropbox
- Laptop can be used with the mouse pad or via touchscreen.

Specifications

- Detector Technology: Indirect conversion, a-Si TFT
- Scintillator: CSI
- Dimensions: W x H x D: 327 mm x 268 mm x 16 mm (12.9" x 10.6" x 0.6")
- Active Area: 233.47 x 291.84 mm (9.2" x 11.5")
- Pixel Matrix: 3,072 x 3,840 pixel
- Resolution (lp/mm): 3.6
- AD Conversion 16 Bit
- Pixel pitch: 76 µm
- Data Transmission: GigE
- Trigger Mode: AED (Automatic Exposure Detection) / external trigger
- Preview Image Time (S) <2
- Full Image Time (S) 5
- Housing Material: Magnesium aluminium alloy
- Weight Limits: Local Load: 100kg at Ø 40 mm / Uniform load: 170kg at effective image area overall
- Battery type: Lithium Ion polymer
- Charging time: 3 hours
- Working time with one battery: 3.5 hours @ 600 images

ADDITIONAL INFORMATION

Animal	Equine
SKU	XR-EXSY000007-PROMO24
Catagory Sub Heading	10 x 12 Wi-Fi Panel - Mobile DR System
System Type	Mobile
Upgrade Panel	No
Wifi	Wifi
A.E.D. Compatible	Yes

Image Transfer Time

<5s

Computer Type

Laptop

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights