

## EXAMION X-R STATIC CLASSIC COMPLETE DR SETUP

XR-EXY000008-PROMO25-BUNDLE\*



### PRODUCT DESCRIPTION

**A complete, all-in-one radiography system ideal for small animal practices wanting to upgrade to DR technology. Simply plug and go...**

Designed to meet all the requirements of a small animal X-ray system. The Examion Static Classic provides you with one complete setup, allowing you to take the image and see the results instantly displayed on the large touchscreen.

The system really is as simple as plug and go, plugging directly into your mains socket, no expensive electrical power connections are needed for installation, and the onboard computer is powered by the press of a button. Bringing all of your controls, settings and analytical tools to your fingertips on the touchscreen monitor.

This complete small animal radiography has been designed to provide superior image quality on an incredibly user-friendly, all-in-one system. For quick and efficient diagnosis.

Powered by Examion AQS software, you can choose patient pre-sets from the main screen, select from the pre-set organ/placement options and simply use the floating table top with foot pedal to position your patient in the right place, all whilst viewing your images on the screen. The software also offers analytical tools such as measurements and image comparisons and allows you to archive the images directly to your practice PACS system.

#### **This package includes as standard:**

- XL (17" x 17") DR Wi-Fi Panel
- Static Classic Unit
- Lenovo PC
- Eizo Multi-touch Monitor

- The following AQS modules are already installed: PACS vet for image storage, 1 x license for AQS.

### **This package also currently includes the following promotional extras:**

- Extra Send/store AQS module - enabling you to send images from your x-ray system directly to a third-party PACS system, such as your practice software for patient records.
- Extra AQS WEB VET module - for online viewing from any device by logging in to your license.
- 5-year panel warranty
- 1st-year AQS software maintenance contract
- 2 x extra AQS licenses for use by multiple users across different locations.

### **Optional Upgrades:**

**Orthopedic (Multi-Room) Surgery Upgrade** - The XL panel installed in the Static Classic system can be used mobile throughout the clinic, with any x-ray tube. This add on supplies you with a second notebook style laptop for easy use away from the Static Classic station. Providing you with not just a complete static small animal solution but also a mobile radiography package.

**Premium Panel Upgrade** - Upgrade to the latest technology, with the premium 100-micron panel, for the highest quality images available.

### **Static Classic Key Features:**

- All-in-one system, designed for Radiography rooms in small animal practices
- All powered by one mains cable - simply plug and go - no extra costly power connections needed for installation.
- Floating top for easy animal positioning powered by a foot pedal.
- Advanced CsI scintillator technology ensures outstanding x-ray images.
- Powered by Examion AQS PACS
- Images displayed in <4 seconds
- Animal-specific organ programs for setting the system and detector to the perfect setup with just one click.
- High power generator for no-compromise quality
- 22" Multi-Touchscreen monitor - zoom in or rotate images by simply spreading your fingers.
- Optional connection to your surgery systems
- Designed to substantially improve your radiography quality and workflow.
- Optional Wi-Fi panel enabling you to remove the panel from the Static Classic and use it around the practice.
- Fully supported by Burtons & Examion - Software issues can be diagnosed and fixed remotely, reducing x-ray equipment downtime.
- AQS software includes pre-sets for easy acquisition, as well as image archive and viewer, with analytical tools.
- Able to add extra Examion bolt-on modules as your practice radiography needs grow.

### **AQS Software:**

- 30 years of continuous software design and development
- Expandability built into the core of the system.
- Intuitive user interface for easy image capture and manipulation
- Highest-quality images from your chosen digital device
- Multi-option search feature to locate 'lost' images
- Future-proof design provides upgrade options to the latest innovations in CR and DR technology
- Professional PACS application with integrated PACS server
- Image capture from a range of CR and DR devices
- Connect to multiple CR/DR devices within the same practice
- Veterinary species-specific processing algorithms
- DICOM 3.0 conformant
- DICOM worklist communication where available in PIMS (PMS) applications
- Can act as a DICOM modality
- Expandable via modules to communicate with, and store images, and 'movie' clips from CT, MRI, Ultrasound and other imaging modalities
- Expandable via modules to receive and store DICOM, JPG, BMP, and other common image files

- Expandable image storage space via a separate PACS archive server
- Multi-option search function
- Multi-monitor viewing support
- Patient CD creation with integral viewing software
- Email images as standard attachments
- Standard annotation functions, linear length, scaled linear length angle, open-angle, area, textbox, short report and more
- Advanced veterinary measuring functions, VHS, Laminitis and more
- Built-in image archiving facility

**Specifications:**

- Generator: 32kw
- Exposure time: 1ms - 2.5s
- Simple wall plug connection
- Table top size : 160cm x 75cm
- Footswitch control for floating table

**DR System:**

- Latest automatic exposure detection flat panel detector
- Rapid image capture and exposure times of 1ms-2.5s
- 17x17in panel size
- Optimised image processing for veterinary application
- 22" Touchscreen monitor

**X-AQS PACS Vet**

- Intuitive X-AQS software
- Import/ export Bitmap & DICOM images
- 351 pre-installed programs
- Server = 1TB
- Reporting in fully functional view screen

**Monitor:**

- 22"
- Colour
- Multi-touch
- Resolution: 1920 x 1080 pixels

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	XR-EXY000008-PROMO25-BUNDLE*
<b>Overall Product Dimensions (mm)</b>	Table height - 78 cm - Table top surface 160 cm x 75 cm (63" x 29.5"), offset ± 35 cm (13.8")
<b>Brand</b>	Examion
<b>Catagory Sub Heading</b>	<b>Static</b> Small Animal Solution
<b>System Type</b>	Station
<b>Upgrade Panel</b>	Yes
<b>Wifi</b>	Full Wifi
<b>A.E.D. Compatible</b>	Yes
<b>Computer Type</b>	Tower

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights