

**BURTONS M880B HANDHELD CAPNOGRAPH AND PULSE OXIMETER  
- CLEARANCE**

129-870-CL-2

**PRODUCT DESCRIPTION**

**Ultra clear, ultra lightweight touchscreen Capnograph & Pulse Oximetry monitor. The perfect addition to any practice!**

Boasting a large 4.3inch, high resolution, auto- rotation touchscreen and a large capacity lithium battery. The M880B is designed for continuous monitoring or spot-checking of CO2 and SPO2 signals in patients.

**Key Features:**

- Extremely light, for carrying and handling Only 0.5KG with full configuration, including battery.
- 4.3inch TFT high resolution with touch screen & auto-rotation for user ease of use with any orientation.
- Protection against ingress of liquid (IPX2 Waterproof)
- Specially designed water filter, ensuring humid environments don't effect readings
- Configurable for flux for horse, dog and cat
- Large capacity lithium battery for up to 18 hours use
- CO2 start/off button to save power
- Selectable layout modes including font and waveforms
- Visual and audio alarm for exceeding limits and Apnea; pitch tone

|                         |                 |
|-------------------------|-----------------|
| Size:                   | 4.3inch         |
| Resolution:             | 480*272         |
| Touch:                  | Resistive Touch |
| Auto-rotation Direction | 4 Direction     |

|                    |  |
|--------------------|--|
| Power Adaptor:     | Input: AC 100 V±2 240V<br>Output: DC 5V/2A         |
| <b>CO2:</b>        |  |
| Measurement Way:   | Infrared Spectrum                                  |
| Measurement Range: | 0-19.7% (0-150mmHg 0-20kPA)                        |
| CO2 Accuracy:      | <5.0% CO2 ±2mmHg<br>>5.0% CO2 <6% of the reading   |
| CO2 Resolution:    | 0.1mmHg  |
| Gas Flow Rate:     | 60-80ml/min  |
| Response Time:     | 1/4±3s   |
| Calibration:       | Offset calibration, auto/ manual; gain calibration |
| <b>SpO2</b>        |  |
| Measurement Range: | 0-100%   |
| Resolution:        | 1%   |
| Accuracy:          | 70-100%± 2%<br>35-69%± 3%<br>0-34% definition      |

**S/N: M67B000009**

## ADDITIONAL INFORMATION

|                           |                                    |
|---------------------------|------------------------------------|
| <b>Animal</b>             | Small Animal                       |
| <b>SKU</b>                | 129-870-cl-2                       |
| <b>Colour/Style</b>       | White                              |
| <b>Brand</b>              | BLT                                |
| <b>Monitor Parameters</b> | CO <sup>2</sup> , SpO <sup>2</sup> |
| <b>Screen Type</b>        | Touchscreen                        |
| <b>Power Source</b>       | Battery                            |
| <b>Optional Modules</b>   | N/A                                |

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights