


HIGH FLOW ARTHROSCOPY SHEATH - CLEARANCE BUNDLE

END-S140-2434R-CL



PRODUCT DESCRIPTION

High Flow Arthroscopy Sheath

1 Rotating Stopcock, 30 Degree For T141-2430 Including Blunt & Sharp Obturators

Includes:

- Wide Angle 2.4mm 30 Degree Arthroscope 14cm WL, Autoclavable (END-T141-2430)
- High Flow Arthroscopy Sheath, 1 Rotating Stopcock, 30 Degree For T141-2430 Including Blunt & Sharp Obturators (END-S140-2434R)

S/N: 0000765228

ADDITIONAL INFORMATION

| | |
|--------|-------------------|
| SKU | END-S140-2434R-cl |
| Animal | Mixed |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights