

## 3.5V PNEUMATIC OTOSCOPE HEAD WITH SPECULA

WA20260-BUNDLE\*



### PRODUCT DESCRIPTION

## 3.5V Welch Allyn Halogen HPX Veterinary Pneumatic Otoscope

The Welch Allyn veterinary pneumatic otoscope provides bright, white light output and a rotating ultraseal lens.

#### Key Features:

- Halogen HPX lamp provides light output for true tissue color and consistent, long-lasting illumination
- Rotating lens with ultraseal enables easy pneumatic otoscopy
- Fiber optic illumination produces cool light

Complete with three polypropylene specula, sizes 4, 5, and 7 mm

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	WA20260-bundle*
<b>Generic Size</b>	X-Small
<b>Colour/Style</b>	Silver
<b>Brand</b>	Welch Allyn
<b>Catagory Sub Heading</b>	Add a handle

Digital

N/A

Diagnostic Head Type

Pneumatic Otoscope

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights