

3.5V CO-AXIAL OPHTHALMOSCOPE WITH LED HEAD

WA11720-LBUNDLE*



PRODUCT DESCRIPTION

Halogen HPX Welch Allyn 3.5 V Coaxial Ophthalmoscope With LED Head

Key Features

- Coaxial optical system for shadow-free spot
- Crossed linear polarizing filter to reduce corneal reflection
- 18 possible aperture combinations for greater versatility
- Red-free filter for use with any aperture
- Long-lasting LED illumination
- Sealed optics for better dust protection

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	WA11720-Lbundle*
Generic Size	X-Small
Colour/Style	Black
Brand	Welch Allyn
Catagory Sub Heading	Add a handle
Digital	N/A
Diagnostic Head Type	Co-axial Ophthalmoscope

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights